

# EmStat MUX8-R2™

POTENTIOSTAT WITH INTEGRATED MULTIPLEXER



## Contents

Description .....	3
Stacking up to 128 channels .....	3
Configurations .....	4
Connectors .....	4
Cell Connections .....	5
Option A (default): .....	5
Option B: .....	5
Option C: .....	5
Option D: .....	5
Supported Switching Modes .....	6
Specifications of General Parameters .....	7
System Specifications .....	8
Potentiostat .....	8
Electrometer .....	8
Integrated MUX8-R2 Multiplexer .....	8
Other .....	8
Functional Diagram .....	9
Dimensions .....	10
PSTrace: Software for Windows .....	11
System requirements .....	11
Also available as extension .....	12

## Description

The EmStatMUX8-R2 potentiostat with multiplexer is designed for up to 8 channels with 2- or 3-electrode sensors or cells. The instrument consists of the MUX8-R2 multiplexer with an integrated EmStat3 or EmStat3+ potentiostat. Specifications can be found on page 8.

## Stacking up to 128 channels

The EmStatMUX8-R2 has a Link connector which can be used to daisy chain to a MUX8-R2 multiplexer, expanding the number of channels. A maximum of 16 multiplexers can be connected in a daisy chain, giving a maximum of 128 channels.

The PStTrace software detects automatically how many multiplexers are daisy chained and shows the available number of channels in the user interface.



**See page 11 for more information about PStTrace.**



*Magnetic feet for easy stacking*

## Configurations

The EmStatMUX8-R2 multiplexer is designed for use up to 128 channels with 2- or 3- electrode sensors or cells.

The multiplexer can be used with different electrode or sensor configurations:


- 1 Eight separate cells or sensors each with a working/sense, reference and counter electrode
- 2 Eight separate cells or sensors each with a working/sense and combined reference and counter electrode
- 3 Cell or sensor array with eight working/sense electrodes sharing one reference and one counter electrode
- 4 Cell or sensor array with eight working/sense electrodes sharing one combined reference/counter electrode

In all configurations the cells can be multiplexed, leaving the non-selected working electrodes either at open circuit (individually floating) or at Ground potential.

When in configurations 3 and 4, the unselected channels are switched to Ground, they will have the working electrode's potential. This is due to the fact that the active WE is always at Ground potential.

You can easily change the hardware configuration of the MUX8-R2 as part of the measurement settings in our PStTrace software or the PStouch app for Android.

*See page 11 for more information about PStTrace.*



*Hardware settings can be changed in the software*

## Connectors

The EmStatMUX8-R2 has the following connectors:

CONNECTOR	FUNCTION
INPUT	Y-cable connects to both potentiostat sensor connector and (digital) AUX
AUX	Can be used to measure auxiliary input like temperature or pH, and to switch or trigger external hardware using two digital control lines that can be set in PStTrace
LINK	Connects to Input of next multiplexer, for daisy-chaining multiple multiplexers.
USB-C	For connecting to PC or Android device
CHANNEL 1-4	Connects to sensor cables 1-4
CHANNEL 5-8	Connects to sensor cables 5-8


*See next page for cable and cell connection options.*

## Cell Connections

### Option A (default):

The channels are divided in two sets of four sensor cables joined with a D-Sub connector.

Order code: [CBL-MUX08R2-SNS-5S]



### Option B:

The cable here shown at the right can be used in case the multiplexer needs to be connected to a fixed setup by means of soldering or screw-terminals.

Order code: [CBL-HD-MUX08R2]



### Option C:

You can also connect one or two screw-terminals directly in the multiplexer.


Order code: [MUX08R2-ST]



### Option D:

The SPE adapter for our MUX8-R2 multiplexer allows you to connect 8x Screen Printed Electrodes (SPE's). The pitch of the SPE connector is 2.54 mm and compatible with the most popular brands of SPE's.

Order code: [MUX08R2-SPE]



## Supported Switching Modes

In *sequential* mode each channel is set before the next measurement starts. In *alternating* mode, the channels are quickly scanned during each interval time giving a virtual-simultaneous measurement across the selected channels.

Voltammetric techniques:	Supported Switching Mode	
	Sequentially	Alternatingly
▪ Linear Sweep Voltammetry	✓	
▪ Cyclic Voltammetry	✓	
▪ Fast Cyclic Voltammetry	✓	
▪ AC Voltammetry	✓	
▪ Differential Pulse Voltammetry	✓	
<b>Pulsed techniques:</b>		
▪ Square Wave Voltammetry	✓	
▪ Normal Pulse Voltammetry	✓	
▪ Stripping Chronopotentiometry	✓	
<b>Amperometric techniques</b>		
▪ Chronoamperometry	✓	✓
▪ Zero Resistance Amperometry	✓	✓
▪ Multistep Amperometry	✓	
▪ Fast Amperometry	✓	
▪ Pulsed Amperometric Detection	✓	
▪ Multiple-Pulse Amperometric Detection	✓	
<b>Galvanostatic techniques</b>		
▪ Linear Sweep Potentiometry	✓	
▪ Chronopotentiometry	✓	✓
▪ Multistep Potentiometry	✓	
▪ Open Circuit Potentiometry	✓	✓
▪ Stripping Chronopotentiometry	✓	
<b>Other</b>		
▪ Mixed Mode	✓	
▪ Impedance Spectroscopy (EIS/GEIS)	✓	

## Specifications of General Parameters

### General pretreatment

Apply conditioning, deposition or begin potential for: 0 – 1600 s

### General voltammetric parameters

Potential range for EmStat3: -3.000 V to +3.000 V

Potential range for EmStat3+: -4.000 V to +4.000 V

Step potential: 0.125 mV to 250 mV

Pulse potential: 0.125 mV to 250 mV

### Limits of some technique specific parameters for EmStat3 and EmStat3+

<b>NPV and DPV:</b>	Scan rate:	0.025 mV/s (0.125 mV step) to 50 mV/s (5 mV step)
	Pulse time:	5 ms to 300 ms
<b>SWV<sup>1</sup>:</b>	Frequency:	1 Hz to 500 Hz <sup>1</sup>
<b>LSV and CV:</b>	Scan rate:	0.01 mV/s (0.1 mV step) to 5 V/s (5 mV step)
	<b>AD:</b>	Interval time: 1 ms to 300 s Run time: 1 s to hours
<b>PAD:</b>	Interval time:	50 ms to 300 s
	Pulse time:	1 ms to 1 s
	Run time:	10 s to hours
<b>MPAD:</b>	Pulse times:	100 ms to 2 s
	Run time:	10 s to hours
	Number of potential levels:	3
<b>Potentiometry at open circuit (OCP):</b>	Interval time:	1 ms to 30 s
	Maximum run time:	hours
<b>Multistep Amperometry:</b>	Interval time:	1 ms to 30 s
	Number of potential levels:	1 to 255
	Number of cycles:	1 to 20000
	Maximum run time:	hours

<sup>1</sup> PSTrace provides the option to measure forward and reverse currents separately.

Note: some limits of parameters are set for practical reasons and can be modified on request.

## System Specifications

	With potentiostat version <b>EmStat3™</b>	With potentiostat version <b>EmStat3+™</b>
▪ dc-potential range	± 3.000 V	± 4.000 V
▪ compliance voltage	± 5 V	± 8 V
▪ applied dc-potential resolution	0.1 mV	0.125 mV
▪ applied potential accuracy	≤ 0.2 % with max. 2 mV offset error	≤ 0.3 % with max. 3 mV offset error
▪ current ranges	1 nA to 10 mA (8 ranges)	1 nA to 100 mA (9 ranges)
▪ maximum measured current	± 20 mA typical and ± 15 mA minimum	± 100 mA typical

### Potentiostat

- current resolution 0.1 % of current range  
1 pA on lowest current range
- current accuracy ≤ 1 % of current range at 1 nA  
≤ 0.5 % at 10 nA  
≤ 0.2 % at 100 nA to 100 uA  
≤ 0.5 % at 1 mA, 10 mA and 100 mA  
all with max. 0.2 % offset error

### Electrometer

- electrometer amplifier input > 100 Gohm // 4 pF
- rise time approx. 100 μs

### Integrated MUX8-R2 Multiplexer

- number of channels 8 (up to 128 channels when daisy chained)
- multiplexer switches 8 x (WE, S, RE and CE)
- on resistance for WE 1.5 ohm typical
- charge injection on WE 20 pC typical
- leakage current < 20 pA (5 pA typical) at 25 °C
- switching time 2 ms


### Other

- housing aluminium: 138 mm x 121 mm x 37 mm
- weight +/- 250 g
- temperature range 0° C to +40° C
- power supply USB
- communication USB-C

*See page 5 for cable and cell connection options*




### Functional Diagram



## Dimensions

Dimensions in mm:



## PSTrace: Software for Windows

Select current ranges for auto ranging and the starting current range.

Switch between plots if curves with different units are available.

The screenshot displays the PSTrace software interface. On the left, the 'Method Editor' panel shows 'Cyclic Voltammometry' settings, including 'Select current range(s)' with options from 100 pA to 100 mA. The main plot area shows a cyclic voltammogram with a peak at 15.530 (1) and 16.774 (1). On the right, the legend shows 'Cyclic Voltammometry [1]\*' and 'Scan 1\*', 'Scan 2\*', and 'Scan 3\*'. Below the plot, there are panels for 'Selected measurement' and 'Curve' properties.

Click on a measurement in the legend to see the available data and to generate more curves.

Click on a curve in the legend to change its title or appearance.

## Other functions in PSTrace 5

- Equivalent Circuit Fitting
- Scripting
- Open your data in Origin and Excel with one click of a button
- Save all available curves, measurement data and methods to a single file
- Browse measurements on PalmSens4's internal storage
- Direct feedback on method parameters

The screenshot shows the 'Equivalent Circuit Analysis - Fitted potential at 100 mV step' window. It features a circuit diagram with components R1, R2, C1, and C2. To the right, there is a plot of 'Imaginary Impedance (Z') / Real Impedance' versus frequency. Below the plot is a table of results:

Element	Fitted Value	Unit	Error%
R 1	198.5	$\Omega$	0.183
R 2	1.456E+4	$\Omega$	0.208
C 1	1.027	$\mu F$	0.038
R 3	9987	$\Omega$	0.221
C 2	0.032	$\mu F$	0.168

### Integration with third party software:

- Excel
- Origin
- Matlab
- ZView



### System requirements

- Minimum PC requirements are:
- Windows 7, 8, or 10 (32-bit or 64-bit)
  - 1 GHz or faster 32-bit (x86) or 64-bit (x64) processor
  - 1 GB RAM (32-bit) or 2 GB RAM (64-bit)

For more information about software visit [www.palmsens.com/software](http://www.palmsens.com/software)



## Also available as extension

The MUX8-R2 multiplexer is also available as generic multiplexer for PalmSens3, PalmSens4 and EmStat Blue.

For more information:

[www.palmsens.com/product/mux8-r2/](http://www.palmsens.com/product/mux8-r2/)

Please don't hesitate to contact PalmSens for more details:  
[info@palmsens.com](mailto:info@palmsens.com)

**PalmSens BV**  
The Netherlands  
[www.palmsens.com](http://www.palmsens.com)

### DISCLAIMER

Changes in specifications and typing errors preserved.  
Every effort has been made to ensure the accuracy of this document. However, no rights can be claimed by the contents of this document.