




Basic **thick film** sensors

The standard dimensions of BST sensors are 25.4 mm length x 7 mm width x 0.635 mm thickness. The standard substrate is aluminium oxide ceramic.

Adaptation of both dimension and substrate is possible upon customer request.



Art.- No.	Article
100 000000	Thick film sensor BST2-WE-C001, RE-Ag/AgCl
120 000000	Thick film sensor BST3-WE-Pt, CE-Pt, RE-Ag/AgCl
121 000000	Thick film sensor BST3-WE-Au, CE-Pt, RE-Ag/AgCl
122 000000	Thick film sensor BST3-WE-C001, CE-Pt, RE-Ag/AgCl
154 000000	Thick film sensor BST2-WE-Au, RE-Ag/AgCl
1410 000000	Flow cell for thick film electrodes, BST2 (two electrode system)
1411 000000	Flow cell for thick film electrodes, BST3 (three electrode system)