




ITALSENS CORROSION CELL SET

Contents

This ItaSens Corrosion Cell Kit will provide you a standard corrosion cell with a sample holder for metal coupons exposing 1 cm diameter round area to the solution as well as the equipment for the Stern-Makrides arrangement. This cell is suitable for many common corrosion related measurements.

#	Quantity	Description
1	2x	Counter Electrode
2	6x	O-ring
3	1x	PTFE ring
4	2x	Threaded Adapter for T24/40 Joint
5	1x	Reference Electrode
6	1x	Reference Electrode Bridge Tube
7	5x	Reference Electrode Bridge Core
8	1x	Purge Tube
9	1x	Electrode Mounting Rod

#	Quantity	Description
10	1x	Screw Nut
11	1x	PTFE Pad
12	1x	Working Electrode Holder
13	5x	Replacement Teflon Gaskets for Mounting Sample
14	1x	Working Electrode Tip
15	1x	1000 ml Corrosion Flask
16	1x	Ball Clip
17	1x	Flat Specimen Holder Kit



Please do not hesitate to contact PalmSens for more details: info@palmSens.com

PalmSens BV
The Netherlands
www.palmSens.com

DISCLAIMER

Changes in specifications and typing errors reserved. Every effort has been made to ensure the accuracy of this document. However, no rights can be claimed by the contents of this document.