


# EmStat4M™ DEVELOPMENT BOARD TOP VIEW



CD-ES4MDEV-11/21-A

# EmStat4M™ DEVELOPMENT BOARD BOTTOM VIEW

## Arduino Headers



Compatible with  
Arduino MKR Series

AREF	+5V
DAC0/A0	VIN
A1	+3V3
A2	GND
A3	RESET
A4	TX
A5	RX
A6	SCL
-0	SDA
-1	MISO
-2	SCK
-3	MOSI
-4	7
-5	6

MKR\_Trig

SW6



MKR to ES

RX  
TX  
D4

## CON10

JST PH (2.0mm)  
Li-Ion Battery Connector

## Battery Power



## RTC Battery Compartment

CR1220 3V  
Real Time Clock  
backup source

## Development Board Manual

A manual for the EmStat4M Development Board is available for more detailed information and instructions. It can be found in the Download section on:  
[www.palmsens.com/es4m-dev](http://www.palmsens.com/es4m-dev)

