



MUX8-R2™

CONNECTOR PIN-OUTS



(male front view)



(female front view)

Channel 1-4

Pin	Function	Pin	Function	Pin	Function
1	CE4	11	RE1	21	AGND
2	RE4	12	RE_SHIELD1	22	SENSE2
3	RE_SHIELD4	13	NC	23	WE1
4	CE3	14	WE4	24	AGND
5	RE3	15	AGND	25	SENSE1
6	RE_SHIELD3	16	SENSE4		
7	CE2	17	WE3		
8	RE2	18	AGND		
9	RE_SHIELD2	19	SENSE3		
10	CE1	20	WE2		

Channel 5-8

Pin	Function	Pin	Function	Pin	Function
1	CE8	11	RE5	21	AGND
2	RE8	12	RE_SHIELD5	22	SENSE6
3	RE_SHIELD8	13	NC	23	WE5
4	CE7	14	WE8	24	AGND
5	RE7	15	AGND	25	SENSE5
6	RE_SHIELD7	16	SENSE8		
7	CE6	17	WE7		
8	RE6	18	AGND		
9	RE_SHIELD6	19	SENSE7		
10	CE5	20	WE6		