

Starter Kit for Electrochemistry

With this set of electrochemical equipment all you need are chemicals and solutions to start your electrochemical experiments.

With the functional, convenient and high quality equipment of BVT and the portable, economical and easy to use research grade potentiostats of PalmSens a bundle to start with electrochemistry right away is now available.

Customers have different requests for their electrochemical experiments. The Starter Kit for Electrochemistry is modular and adjustable to the needs of different customers.



EmStat Blue
potentiostat

OR



PalmSens3
potentiostat/galvanostat/EIS



Controlled stirrer

+




Glass cell with stoppers



Classical electrodes

OR



Screen printed electrodes



Contents of Electrochemistry Kit

Item	Order code	In classic electrode kit	In sensor kit
PalmSens or EmStat Blue	PS3EIS or ESB3	•	•
T4 cell	BVT-TC4	•	•
3 stoppers	Stopper	•	•
stirrer ST1		•	•
BVT Sensor connector	BVT-KA1C.3		•
BVT Starter Kit			•
BVT classic WE	BVT-WCE(S or 1)D(1-4).ES	•	
BVT classic RE	BVT-RCEc.R1.E3	•	
BVT classic CE	BVT-ACEc.ES	•	

For more information visit
www.bvt.cz or www.palmsens.com